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
- **Hepatitis-B virus and HIV**

Results:

- a) In presence of low organic load: 1 minute
- b) In presence of high organic load: 1 minute

Certificate: Prof. Dr. G. Frösner, München, 16 August 1987 / 20 May 1996

- **Rota virus**

Result according to BGVV and DVV guidelines

- a) Without soil load: 1 minute

Certificate: Dr. J. Steinmann, Bremen, 24 November 1994

- **Vaccinia virus**

Result according to BGA (now RKI) and DVV guidelines

- a) Without soil load: 1 minute
- b) With soil load: 1 minute

Certificate: Dr. J. Steinmann, Bremen, 25 August 2005

- **Tested against BVDV (surrogate of Hepatitis-C virus (HCV))**

Result according to BGA (now RKI) and DVV guidelines

- a) Without soil load: 1 minute

Certificate: Dr. J. Steinmann, Bremen, 10 September 2005

Based on the BVDV and Vaccinia virus expertise's (according to RKI / DVV) the product can be claimed as limited virucide (effective against enveloped viruses) according to Robert Koch Institute (RKI), Germany.

▪ Bakterizide / Fungizide Wirkung

Tested according to DGHM (German Society for Hygiene and Microbiology), now VAH

- **Surface Test**

Disinfection of hard surfaces in hospital and general practice according DGHM standard methods (1. September 2001) and catalog of requirements (4. February 2002):

Results:

- a) In presence of low organic load: 1 minute
- b) In presence of high organic load: 1 minute

Certificates: Prof. Dr. R. Schubert, Frankfurt (M), 5 May 2003 / 29 November 2004

- **Surface Test, DGHM VIIth list and ÖGHM guidelines**

Disinfection of hard surfaces in hospital and general practice:

Results: 1 minute

Certificate: Prof. Dr. R. Schubert, Frankfurt (M), 10 February 1995 and 3 September 1999

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- **Surface Test, DGHM VIIth list**

Disinfection of hard surfaces in hospital and general practice:

Result: 1 minute

Certificate: Prof. Dr. E. Marth, Graz, 27 March 1995 and 22 September 1999

- **Modified DGHM guidelines**

Disinfection of acrylic glass or similar surfaces:

Result: 1 minute

Certificate: Prof. Dr. R. Schubert, Frankfurt (M)

- **Qualitative suspension test, DGHM guidelines**

Test strains: Staphylococcus aureus, Escherichia coli, Proteus mirabilis, Pseudomonas aeruginosa, Candida albicans

Result: 5 minutes

Certificate: Prof. Dr. R. Schubert, Frankfurt (M)

- **Untreated Wood, DGHM guidelines**

Disinfection on untreated wood, contaminated with fungi:

Result: 1 hour

Certificate: Prof. Dr. R. Schubert, Frankfurt (M)

Test according to DVG (German Veterinary Medical Society)

- **Food sector**

Results:

a) In presence of low organic load: 30 minutes

b) In presence of high organic load: 30 minutes

Certificates: Prof. Dr. G. Reuter, Berlin, 23 January 1987

Prof. Dr. R. Schubert, Frankfurt (M), 22 July 1992

Test according to European Norms (EN)

- **EN 1276**

Bactericidal results in presence of high / low organic load (Albumin)

Test strains: P. aeruginosa, S. aureus, E. coli and E. hirae

Results:

a) 3,0 g/l Albumin 5 minutes

b) 0,3 g/l Albumin 30 seconds

Certificates: LONZA Basel, Laboratory LPC, 23 May 2000 / 28 May 2004 / 11 April 2005

- **EN 1040**

Bactericidal test, strains: P. aeruginosa, S. aureus

Result: 5 minutes

Certificate: LONZA Basel, Laboratory LPC, 12 May 2000

- **EN 13697 (Surface test)**

Bactericidal result in presence of organic load (Albumin)

Test strains: P. aeruginosa, S. aureus, E. coli and E. hirae

Result:

a) 0,3 g/l Albumin 5 minutes

Certificate: LONZA Basel, Laboratory OPC-E, 22 August 2005

Other investigations:

Investigation on formation of cracks on acrylic glass.

Use area: e.g. Health Clubs

Result: Product does not form cracks on acrylic glass.

Certificate: Röhm GmbH, Darmstadt

Use areas:

Hospital – Institutions – Food Industry – Health Clubs

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